


Product Name: Recombinant Human UBE2W Protein
Catalog #: PEH6528



Summary

Name	UBE2W_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96B02
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MASMQKRLQK ELLALQNDPP PGMTLNEKSV QNSITQWIVD MEGAPGTYE GEKFQLLFKF SSRYPFDSPQ VMFTGENIPV HPHVYSNGHI CLSILTEDWS PALSVQSVCL SIISMLSSCK EKRRPPDNSF YVRTCNKNPK KTKWWYHDDT C

image



Background

Alternative Names	Ubiquitin-conjugating enzyme E2 W (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme W) (N-terminal E2 ubiquitin-conjugating enzyme) (EC 2.3.2.25) (N-terminus-conjugating E2) (Ubiquitin carrier protein W) (Ubiquitin-conjugating enzyme 16) (UBC-16) (Ubiquitin-protein ligase W)
--------------------------	--

Note

Product Name: Recombinant Human UBE2W Protein
Catalog #: PEH6528



For Research Use Only , Not for Diagnostic Use