


Product Name: Recombinant Human ABCC11 Protein
Catalog #: PEH6537



Summary

Name	ABCCB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96J66
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FKKFFLQESP VFYVQTLQDP SKALVFEEAT LSWQQTCPGI VNGALELERN GHASEGMTRP RDALGPEEEG NSLGPELHKI NLVSKGMML GVCGNTGSGK SSLLSAILEE MHLLEGSVGV QGSLAYVPQQ AWIVSGNIRE NILMGGAYDK ARYLQVLHCC SLNRDLELLP FGDMTEIGER GLNLSGGQKQ RISLARAVYS DRQIYLLDDP LSAVDAHVGK HIFEECIKKT LRGKTVVLVT HQLQYLEFCG QIILENGKI CENGTHSELM QKKGKYAQLI QKMHKEATSD MLQDTAKIAE KPKVESQALA TSLEESLNGN AVPEHQLTQE EEMEEGSLSW RVYHHYIQAA

image



Background

Alternative Names ATP-binding cassette sub-family C member 11 (Multidrug resistance-associated protein 8)

Product Name: Recombinant Human ABCC11 Protein
Catalog #: PEH6537



Note

For Research Use Only , Not for Diagnostic Use