


**Product Name: Recombinant Human RAD1 Protein**  
**Catalog #: PEH6543**



## Summary

<b>Name</b>	RAD1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-SUMO & N-His
<b>Accession #</b>	O60671
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	44 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	EDDQYSLVAS LDNVRNLSTI LKAIHFREHA TCFATKNGIK VTVENAKCVQ ANAFIQAGIF QEFKVQEEVS TFRINLTVLL DCLSIGSSP MPGTLTALRM CYQGYGYPLM LFLEEGGVVT VCKINTQEPE ETLDFDFCST NVINKIILQS EGLREAFSEL DMTSEVLQIT MSPDKPYFRL STFGNAGSSH LDYPKSDSLM EAFHCNQTVV NRYKISLLKP STKALVLSCK VSIRTDNRGF LSLQYMIRNE DGQICFVEYY CCPDEEVPE

## image



## Background

**Alternative Names** Cell cycle checkpoint protein RAD1 (hRAD1) (EC 3.1.11.2) (DNA repair exonuclease rad1 homolog) (Rad1-like DNA damage checkpoint protein)

**Product Name: Recombinant Human RAD1 Protein**  
**Catalog #: PEH6543**

---



**Note**

For Research Use Only , Not for Diagnostic Use