


Product Name: Recombinant Human R3HDM2 Protein
Catalog #: PEH6548



Summary

Name	R3HD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2K5
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEIMKESEKK LVEESVNKNK FISKTPSKEE IEKECEDTSL RQETQRRTSN HGHARKRAKS NSKLKLVRS L AVCEESSTPF ADGPLETQDI IQLHISCPSD

image



Background

Alternative Names R3H domain-containing protein 2

Note

For Research Use Only , Not for Diagnostic Use