



Product Name: Recombinant Human RTN4RL2 Protein
Catalog #: PEH6549



Summary

Name	R4RL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q86UN3
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLWAWFQRAR VSSSDVTCAT PPERQGRDLR ALREADFQAC PPAAPTRPGS RARGNSSSNH LYGVAEAGAP PADPSTLYRD LPAEDSRGRQ GGDAPTEDDY WGGYGGEDQR

image



Background

Alternative Names	Reticulon-4 receptor-like 2 (Nogo receptor-like 3) (Nogo-66 receptor homolog 1) (Nogo-66 receptor-related protein 2) (NgR2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use