


Product Name: Recombinant Human RGS8 Protein
Catalog #: PEH6550



Summary

Name	RGS8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P57771
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAALLMPRRN KGMRTRLGCL SHKSDSCSDF TAILPDKPNR ALKRLSTEEA TRWADSFVDL LSHKYGVAAF RAFLKTEFSE ENLEFWLACE EFKKTRSTAK LVSKAHRIFE EFVDVQAPRE VNIDFQTREA TRKNLQEPSL TCFDQAQGVK HSLMEKDSYP RFLRSKMYLD LLSQSQRRLS

image



Background

Alternative Names Regulator of G-protein signaling 8 (RGS8)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RGS8 Protein
Catalog #: PEH6550

