


Product Name: Recombinant Human ADGRF2P Protein
Catalog #: PEH6554



Summary

Name	AGRF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8IZF7
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GWVQLDQLPS PNAVSSDQVH CSAGCTHRKC GWAASKSKEK VPARPHGVCD GVCTDYSQCT QPCPPDTQGN MGFSCRQKTW HKITDTCQTL NALNIFEEDS RLVQPFEDNI KISVYTGKSE TITDMLLQKC PTDLSCVIRN

image



Background

Alternative Names Adhesion G-protein coupled receptor F2 (G-protein coupled receptor 111)
(G-protein coupled receptor PGR20)

Note

For Research Use Only , Not for Diagnostic Use