



Product Name: Recombinant Human ARSH Protein
Catalog #: PEH6558



Summary

Name	ARSH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5FYA8
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IIQQPMKEEK VASLMLKEAL AFIERYKREP FLLFFSFLHV HTPLISKKKF VGRSKYGRYG DNVEEMDWMV GKILDALDQE RLANHTLVYF TSDNGGHLEP LDGAVQLGGW NGIYKGGKGM GGWEGGIRVP GIFRWPSVLE AGRVINEPTS LMDIYPTLSY IGGGILSQDR VIDGQNL MPL

image



Background

Alternative Names Arylsulfatase H (ASH) (EC 3.1.6.-)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ARSH Protein
Catalog #: PEH6558

