


Product Name: Recombinant Human ATL1 Protein
Catalog #: PEH6559



Summary

Name	ATLA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WXF7
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KAGPVQVLIV KDDHSFELDE TALNRILLSE AVRDKVAV SVAGAFRKGK SFLMDFMLRY MYNQESVDWV GDYNEPLTGF SWRGGSERET TGIQIWSEIF LINKPDGKKV AVLLMDTQGT FDSQSTLRDS ATVFALSTMI SSIQVYNLSQ NVQEDDLQHL QLFTEYGRLA MEETFLKPFQ SLIFLVRDWS FPYEFSYGAD GGAKFLEKRL KVSGNQHEEL QNVRKHIHSC FTNISCFLLP HPGLKVATNP NFDGKLKEID DEFIKNLIL IPWLLSPESL DIKEINGNKI TCRGLVEYFK AYIKIYQGEE LPHPKSMLQA TAEANNLAAV ATAKDTYNKK MEEICGGDKP FLAPNDLQTK HLQLKEESVK LFRGVKKMGG EEFSTRYLQQ LESEIDELYI

image



Background

Alternative Names Atlastin-1 (EC 3.6.5.-) (Brain-specific GTP-binding protein) (GTP-binding

Product Name: Recombinant Human ATL1 Protein
Catalog #: PEH6559



protein 3) (GBP-3) (hGBP3) (Guanine nucleotide-binding protein 3) (Spastic paraplegia 3 protein A)

Note

For Research Use Only , Not for Diagnostic Use