



Product Name: Recombinant Human ATG9B Protein
Catalog #: PEH6577



Summary

Name	ATG9B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q674R7
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRRLGRWGDL GPGSVPLLPM PLPPPPPPSC RGPGGGRISI FSLSPAPHTR SSPSSFSPPT AGPPCSVLQG TGASQSCHSA LPATPPTQ AQPAMTPASA SPSWGSHSTP PLAPATPTPS QQCPQDSPGL RVGPLIPEQD YERLEDCDPE GSQDSPIHGE EQQPLLHVPE GLRGSWHHIQ NLDSFFTKIY

image



Background

Alternative Names Autophagy-related protein 9B (APG9-like 2) (Nitric oxide synthase 3-overlapping antisense gene protein) (Protein sONE)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ATG9B Protein
Catalog #: PEH6577

