


Product Name: Recombinant Human ABHD16A Protein
Catalog #: PEH6579



Summary

Name	ABHGA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95870
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ENKRQLANYN FDFRSWPVDF HWEEPSSRKE SRGGPSRRGV ALLRPEPLHR GTADTLLNRV KKLPCQITSY LVAHTLGRRM LYPGSVYLLQ KALMPVLLQG QARLVEECNG RRAKLLACDG NEIDTMFVDR RGTAEPQGQK LVICCEGNAG FYEVGCVSTP LEAGYSVLGW NHPGFAGSTG VPFPQNEANA MDVVVQFAIH RLGFQPQDII IYAWSIGGFT ATWAAMSYPD VSAMILDASF DDLVPLALKV MPDSWRGLVT RTVRQHLNLN

image



Background

Alternative Names	Phosphatidylserine lipase ABHD16A (EC 3.1.-.-) (Alpha/beta hydrolase domain-containing protein 16A) (Abhydrolase domain-containing protein 16A) (HLA-B-associated transcript 5) (hBAT5) (Monoacylglycerol lipase
--------------------------	--

Product Name: Recombinant Human ABHD16A Protein
Catalog #: PEH6579



ABHD16A) (EC 3.1.1.23) (Protein G5)

Note

For Research Use Only , Not for Diagnostic Use