


Product Name: Recombinant Human ABCA6 Protein
Catalog #: PEH6580



Summary

Name	ABCA6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N139
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IEKEIDAEHP SDDYFEPVAP EFQGKEAIRI RNVKKEYK GK SGKVEALKGL LFDIYEGQIT AILGHSGAGK SLLNILNGL SVPTEGSVTI YNKNLSEMQR LEEIRKITGV CPQFNVQFDI LTVKENLSLF AKIKGIHLKE VEQEVQRILL ELDMQNIQDN LAKHLSEGQK RKLTFGITIL GDPQILLDE PTTGLDPFSR DQVWSLLRER RADHVILFST QSMDEADILA DRKVIMSNR LKACAGSSMFL KRRWGLGYHL SLHRNEICNP EQITSFITH IPDAKLTEN KEKLVYTLPL ERTNTFPDLF SLDKCS DQG VTGYDISMST LNEVFMKLEG QSTIEQDFEQ VEMIRDSESL NEMELAHSSF

image



Background

Alternative Names ATP-binding cassette sub-family A member 6

Product Name: Recombinant Human ABCA6 Protein
Catalog #: PEH6580



Note

For Research Use Only , Not for Diagnostic Use