



Product Name: Recombinant Human ACSL4 Protein
Catalog #: PEH6584



Summary

Name	ACSL4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60488
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DWQEGGYTIN DKNPRGEIV IGGQNISMGY FKNEEKTAED YSVDENGQRW FCTGDIGEFH PDGCLQIIDR KKDLVKLQAG EYVSLGKVEA ALKNCPLIDN ICAFKSDQS YVISFVVPNQ KRLTLAQKQK GVEGTWVDIC NNPAMEAEIL KEIREAANAM KLERFEIPIK VRLSPEPWTP ETGLVTDAFK LKRKELRNHY LKDIERMYGG

image



Background

Alternative Names	Long-chain-fatty-acid--CoA ligase 4 (EC 6.2.1.3) (Arachidonate--CoA ligase) (EC 6.2.1.15) (Long-chain acyl-CoA synthetase 4) (LACS 4)
--------------------------	---

Note

Product Name: Recombinant Human ACSL4 Protein
Catalog #: PEH6584



For Research Use Only , Not for Diagnostic Use