


Product Name: Recombinant Human ASCL4 Protein
Catalog #: PEH6587



Summary

Name	ASCL4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6XD76
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAFLRKRNER ERQVRVCVNE GYARLRDHLP RELADKRLSK VETLRAAIDY IKHLQELLER QAWGLEGAAG AVPQRRRAECN SDGESKASSA PSPSSEPEEG

image



Background

Alternative Names	Achaete-scute homolog 4 (ASH-4) (hASH4) (Achaete-scute-like protein 4) (Class A basic helix-loop-helix protein 44) (bHLHa44)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use