


Product Name: Recombinant Human EPHA8 Protein
Catalog #: PEH6593



Summary

Name	EPHA8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P29322
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQVVVIRQER AGQTSVSLW QEPEQPNGII LEYEIKYYEK DKEMQSYSTL KAVTTRATVS GLKPGTRYVF QVRARTSAGC GRFSQAMEVE TGKPRPRYDT

image



Background

Alternative Names	Ephrin type-A receptor 8 (EC 2.7.10.1) (EPH- and ELK-related kinase) (EPH-like kinase 3) (EK3) (hEK3) (Tyrosine-protein kinase receptor EEK)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use