


Product Name: Recombinant Human HOXB9 Protein
Catalog #: PEH6601



Summary

Name	HXB9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17482
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YVDSIISHES EDAPPAKFPS GQYASSRQPG HAEHLEFPSC SFQPKAPVFG ASWAPLSPHA SGSLPSVYHP YIQPQGVPPA ESRYLRTWLE PAPRGEAAPG QGQAAVKAEP LLGAPGELLK QGTPEYSLET SAGREAVLSN QRPYGDNKI CEGSEDKERP DQTNPSANWL HARSSRKKRC PYTKYQTLLEL EKEFLFNMYL TRDRRHEVAR LLNLSERQVK

image



Background

Alternative Names	Homeobox protein Hox-B9 (Homeobox protein Hox-2.5) (Homeobox protein Hox-2E)
--------------------------	------------------------------------------------------------------------------

Note

Product Name: Recombinant Human HOXB9 Protein
Catalog #: PEH6601



For Research Use Only , Not for Diagnostic Use