



Product Name: Recombinant Human HOXC9 Protein
Catalog #: PEH6604



Summary

Name	HXC9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P31274
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YVDSLISHDN EDLLASRFPA TGAHPAAARP SGLVPDCSDF PSCSFAPKPA VFSTSWAPVP SQSSVYHPY GPQPHLGADT RYMRTWLEPL SGAVSFPSFP AGGRHYALKP DAYPGRRADC GPGEGRSYPD YMYGSPGELR DRAPQTLPSF EADALAGSKH KEEKADLDPS NPVANWIHAR STRKKRCPYT KYQTLELEKE FLFNMYLTRD RRYEVARVLN LTERQVKIWF QNRRMKMKKM

image



Background

Alternative Names Homeobox protein Hox-C9 (Homeobox protein Hox-3B)

Note

Product Name: Recombinant Human HOXC9 Protein
Catalog #: PEH6604



For Research Use Only , Not for Diagnostic Use