


Product Name: Recombinant Human ICOSLG Protein
Catalog #: PEH6616



Summary

Name	ICOSL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O75144
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QEKEVRAMVG SDVELSCACP EGSRFDLNDV YVYWQTSESK TVVTYHIPQN SSENVDSRY RNRALMSPAG MLRGDFSLRL FNVTPQDEQK FHCLVLSQSL GFQEVLSVEV TLHVAANFSV PVVSAPHSPS QDELTFTCTS INGYPRPNVY WINKTDNSLL DQALQNDTVF LNMRLGYDVV SVLRIARTPS VNIGCCIENTV LLQQNLTVGS QTGNDIGERD KITENPVSTG

image



Background

Alternative Names	ICOS ligand (B7 homolog 2) (B7-H2) (B7-like protein GI50) (B7-related protein 1) (B7RP-1) (CD antigen CD275)
--------------------------	--

Note

Product Name: Recombinant Human ICOSLG Protein
Catalog #: PEH6616



For Research Use Only , Not for Diagnostic Use