


Product Name: Recombinant Human IFNL4 Protein
Catalog #: PEH6620



Summary

Name	IFNL4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	K9M1U5
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAPRRCLLSH YRSLEPRTLA AAKALRDRIE EEALSWGQRN CSFRPRRDPP RPSSCARLRH VARGIADAQA VLSGLHRSEL LPGAGPILEL LAAAGRDVAA CLELARGSS RKVPGAQKRR HKPRRADSPR CRKASVFNL LRLLTWELRL AAHSGPCL

image



Background

Alternative Names Interferon lambda-4 (IFN-lambda-4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IFNL4 Protein
Catalog #: PEH6620

