


Product Name: Recombinant Human NLRP7 Protein
Catalog #: PEH6621



Summary

Name	NALP7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WX94
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QTLLEQLNED ELKSFKSLW AFPLEDLQK TPWSEVEEAD GKKLAEILVN TSENWIRNA TVNILEEMNL TELCKMAKAE MMEDGQVQEI DNPELGDAEE DSELAKPGEK

image



Background

Alternative Names	NACHT, LRR and PYD domains-containing protein 7 (Nucleotide-binding oligomerization domain protein 12) (PYRIN-containing APAF1-like protein 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use