


Product Name: Recombinant Human MYG1 Protein
Catalog #: PEH6624



Summary

Name	MYG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9HB07
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YRDAEIVRTR DPEKLASCDI VVDVGGEYDP RRHRYDHHQR SFTETMSSLS PGKPWQTKLS SAGLIYLHFG HKLLAQLLGT SEEDSMVGTL YDKMYENFVE EVDAVDNGIS QWAEGEPRYA LTTTLSARVA RLNPTWNHPD QDTEAGFKRA MDLVQEEFLQ RLDFYQHSWL PARALVEEAL AQRFAQVDPSG EIVELAKGAC PWKEHLYHLE SGLSPPVAIF FVIYTDQAGQ WRIQCVPKEP HSFQSRPLPL EPWRGLRDEA

image



Background

Alternative Names MYG1 exonuclease (EC 3.1.-.-)

Product Name: Recombinant Human MYG1 Protein
Catalog #: PEH6624



Note

For Research Use Only , Not for Diagnostic Use