



## Summary

<b>Name</b>	NEIL3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8TAT5
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	17 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	TLERKTKQNQ ILDEEFQNSP PASVCLNDIQ HPSKKTNDI TQPSSKVNIS PTISSESKLF SPAHKKPKTA QYSSPELKSC NPGYSNSELQ INMTDGPRTL NPDSPRCSKH NRLCILRVVG KDGENKGRQF YACPLPREAQ

## image



## Background

<b>Alternative Names</b>	Endonuclease 8-like 3 (EC 3.2.2.-) (EC 4.2.99.18) (DNA glycosylase FPG2) (DNA glycosylase/AP lyase Neil3) (Endonuclease VIII-like 3) (Nei-like protein 3)
--------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human NEIL3 Protein**  
**Catalog #: PEH6626**

---

