


Product Name: Recombinant Human NAXE Protein
Catalog #: PEH6627



Summary

Name	NNRE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NCW5
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELFNEYQFSV DQLMELAGLS CATAIAKAYP PTSMSPPT VLVICGPGNN GGDGLVCARH LKLFGEPTI YYPKRPNKPL FTALVTQCQK MDIPFLGEMP AEPMTIDELY ELVVDAIFGF SFKGDVREPF HSILSVLKGL TVPIASIDIP SGWDVEKGNA GGIQPDLLIS LTAPKKSATQ FTGRYHYLGG RFVPPALEKK YQLNLPPYPD

image



Background

Alternative Names	NAD(P)H-hydrate epimerase (EC 5.1.99.6) (Apolipoprotein A-I-binding protein) (AI-BP) (NAD(P)HX epimerase) (YjeF N-terminal domain-containing protein 1) (YjeF_N1)
--------------------------	---

Product Name: Recombinant Human NAXE Protein
Catalog #: PEH6627



Note

For Research Use Only , Not for Diagnostic Use