


Product Name: Recombinant Human NFXL1 Protein
Catalog #: PEH6628



Summary

Name	NFXL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZNB6
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLSCCKNQCP KELPCGHRCK EMCHPGECPF NCNQKVKLRC PCKRIKKELQ CNKVRENQVS IECDDTCKEM KRKASEIKEA EAKAALEEEK RRQQAELEAF

image



Background

Alternative Names NF-X1-type zinc finger protein NFXL1 (Ovarian zinc finger protein) (hOZFP)

Note

For Research Use Only , Not for Diagnostic Use