


Product Name: Recombinant Human LRFN1 Protein
Catalog #: PEH6631



Summary

Name	LRFN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9P244
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLHKLPPDGL FLRSQGTGPK PPTPLTVSFG GNPLHCNCEL LWLRRRLTRED DLETCATPEH LTDRYFWSIP EEEFLCEPPL ITRQAGGRAL VVEGQAVSLR CRAVGDPPEV VHWVAPDGRL LGNSSRTRVR GDGTLDTVIT TLRDSGTFTC IASNAAGEAT APVEVCVVPL PLMAPPPAAP PPLTEPGSSD IATPGRPGAN

image



Background

Alternative Names	Leucine-rich repeat and fibronectin type III domain-containing protein 1 (Synaptic adhesion-like molecule 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LRFN1 Protein
Catalog #: PEH6631

