


Product Name: Recombinant Human NSG2 Protein
Catalog #: PEH6633



Summary

Name	NSG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y328
Host	E.coli
Species	Human
Molecular Mass	8 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGTKPPSVED GFQTVPLITP LEVNHLQLPA PEKVIVKTRT EYQPEQKNKG

image



Background

Alternative Names	Neuronal vesicle trafficking-associated protein 2 (Neuron-specific protein family member 2) (Protein p19) (Hmp19)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use