


Product Name: Recombinant Human TRIT1 Protein
Catalog #: PEH6636



Summary

Name	MOD5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H3H1
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IQGHKPTATP IKMPYNEAEN KRSYHLCDLC DRIIIGDREW AAHIKSKSHL NQLKKRRRLD SDAVNTIESQ SVSPDHNKEP KEGSPGQND

image



Background

Alternative Names tRNA dimethylallyltransferase (EC 2.5.1.75) (Isopentenyl-diphosphate:tRNA isopentenyltransferase) (IPP transferase) (IPPT) (hGRO1) (tRNA isopentenyltransferase 1) (IPTase)

Note

For Research Use Only , Not for Diagnostic Use