


Product Name: Recombinant Human KCNT1 Protein
Catalog #: PEH6637



Summary

Name	KCNT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5JUK3
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQIPLWSQRV IYLGSAKLD QDLMRAKMDN GEACFILSSR NEVDRTAADH QTILRAWAVK DFAPNCPLYV QILKPENKFH VKFADHVVC EECKYAMLAL NCICPATSTL ITLLVHTSRG QEGQESPEQW QRMYGRCSGN EVYHIRMGDS KFFREYEGKS FTYAAFHAHK KYGVCLIGLK REDNKSILLN PGPRHILAAS DTCFYINITK EENSAFIFKQ EEKRKKRAFS GQGLHEGPAR LPVHSIIASM GTVAMDLQGT EHRPTQSGGG GGGSKLALPT ENGSGSRRPS IAPVLELADS SALLPCDLLS DQSEDEVTPS DDEGLSVVEY VKGYPPNSPY IGSSPTLCHL LPVKAPFCCL RLDKGCKHNS YEDAKAYGFK NKLIIVSAET AGNGLYNFIV PLRAYRSRK ELNPVILLDD

image



Background

Product Name: Recombinant Human KCNT1 Protein
Catalog #: PEH6637



Alternative Names

Potassium channel subfamily T member 1 (KCa4.1)

Note

For Research Use Only , Not for Diagnostic Use