


Product Name: Recombinant Human NAALAD2 Protein
Catalog #: PEH6639



Summary

Name	NALD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y3Q0
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WDAEEFGLLG STEWAEENVK ILQERSIAYI NSDSSIEGNY TLRVDCTPLL YQLVYKLTKE IPSPDDGFES KSLYESWLEK

image



Background

Alternative Names	N-acetylated-alpha-linked acidic dipeptidase 2 (EC 3.4.17.21) (Glutamate carboxypeptidase III) (GCPIII) (N-acetylated-alpha-linked acidic dipeptidase II) (NAALADase II)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use