


Product Name: Recombinant Human CA8 Protein
Catalog #: PEH6646



Summary

Name	CAH8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P35219
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YEEGVWGLV FPDANGEYQS PINLNSREAR YDPSLLDVRL SPNYVVCRCDC EVTNDGHTIQ VILKSKSVLS GGPLPQGHEF ELYEVRFWHG RENQRGSEHT VNFKAFPMEL HLIHWNSTLF GSIDEAVGKP HGIAIIALFV QIGKEHVGLK AVTEILQDIQ YKGKSKTIPC FNPNTLLPDP LLRDYWVYEG SLTIPPCSEG VTWILFRYPL TISQLQIEEF RRLRTHVKGA ELVEGCDGIL GDNFRPTQPL

image



Background

Alternative Names	Carbonic anhydrase-related protein (CARP) (Carbonic anhydrase VIII) (CA-VIII)
--------------------------	---

Note

Product Name: Recombinant Human CA8 Protein
Catalog #: PEH6646



For Research Use Only , Not for Diagnostic Use