


Product Name: Recombinant Human GPC2 Protein
Catalog #: PEH6647



Summary

Name	GPC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8N158
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SCAETRQVLG ARGYSLNLIP PALISGEHLR VCPQEYTCSS SETEQRLLIRE TEATFRGLVE DSGSFLVHTL AARHRKFDEF FLEMLSVAQH SLTQLFSHSY GRLYAQHALI FNGLFSRLRD FYGESGEGLD DTLADFWAQL LERVFPLLHP QYSFPPDYLL CLSRLASSTD GSLQPFGDSP

image



Background

Alternative Names Glypican-2 [Cleaved into: Secreted glypican-2]

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GPC2 Protein
Catalog #: PEH6647

