


Product Name: Recombinant Human Ig Kappa Light Chain Protein
Catalog #: PEH6658

Summary

Name	IGK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P0DOX7
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DIQMTQSPST LSASVGDRVT ITCRASQSIN TWLAWYQQKP GKAPKLLMYK ASSLESGVPS RFIGSGSGTE FTLTISLQP DDFATYYCQQ YNSDSKMFGQ GTKVEVKGTV AAPSVFIFPP SDEQLKSGTA SVVCLLNNFY PREAKVQWKV DNALQSGNSQ ESVTEQDSKD STYLSSTLT LSKADYEKHK VYACEVTHQG LSPVTKSFN RGEK

image



Background

Alternative Names Immunoglobulin kappa light chain (Immunoglobulin kappa light chain EU)

Note

Product Name: Recombinant Human Ig Kappa Light Chain Protein
Catalog #: PEH6658



For Research Use Only , Not for Diagnostic Use