


Product Name: Recombinant Human MED31 Protein
Catalog #: PEH6662



Summary

Name	MED31_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y3C7
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DAGNRLRFQL ELEFVQCLAN PNYLNFLAQR GYFKDKAFVN YLKYLWYWKD PEYAKYLYKYP QCLHMLELLQ YEHFRKELVN AQCAKFIDEQ QILHWQHYSR KRMRLQQALA EQQQQNNTSG

image



Background

Alternative Names Mediator of RNA polymerase II transcription subunit 31 (Mediator complex subunit 31) (Mediator complex subunit SOH1) (hSOH1)

Note

For Research Use Only , Not for Diagnostic Use