


Product Name: Recombinant Human DLGAP3 Protein
Catalog #: PEH6663



Summary

Name	DLGP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95886
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAEGPARRCS SADGLDGPAM GARTLELAPV PPRASPKPPT LIIKTIPGRE ELRSLARQRK WRPSIGVQVE TISDSDTENR SRREFHSIGV QVEEDKRRAR

image



Background

Alternative Names	Disks large-associated protein 3 (DAP-3) (PSD-95/SAP90-binding protein 3) (SAP90/PSD-95-associated protein 3) (SAPAP3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use