


Product Name: Recombinant Human CYB5B Protein
Catalog #: PEH6676



Summary

Name	CYB5B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43169
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ETSVTYRLE EVAKRNSLKE LWLVHGRVY DVTRFLNEHP GGEEVLLEQA GVDASESFED VGHSSDAREM LKQYYIGDIH PSDLKPEGS KDPSKNDTCK

image



Background

Alternative Names	Cytochrome b5 type B (Cytochrome b5 outer mitochondrial membrane isoform)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use