


Product Name: Recombinant Human KEL Protein
Catalog #: PEH6677



Summary

Name	KELL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P23276
Host	E.coli
Species	Human
Molecular Mass	87 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VCLDLRDHYL ASGNTSVAPC TDFSFCACGR AKETNNSFQE LATKNKNRLR RILEVQNSWH PGSGEEKAFQ FYNSCMDTLA IEAAGTGPLR QVIEELGGWR ISGKWTSNLF NRTLRLMSQ YGHFPFFRAY LGPHPASPHT PVIQIDQPEF DVPLKQDQEQ KIYAQIFREY LTYLNQLGTL LGGDPSKVQE HSSLISITS RLFQFLRPLE QRRAQGKLFQ MVTIDQLKEM APAIDWLSCL QATFTPMSLS PSQSLVVHDV EYLNMSQLV EEMLLKQRDF LQSHMILGLV VTLSPALDSQ FQEARRKLSQ KLRELTEQPP MPARPRWMKC VEETGTFFEP TLAALFVREA FGPSTRSAAM KLFTAIRDAL ITRLRNLPWM NEETQNMAQD KVAQLQVEMG ASEWALKPEL ARQEYNDIQL GSSFLQSVLS CVRSLRARIV QSFLQPHPQH RWKVSPWDVN AYYSVSDHVV VFPAGLLQPP FFHPGYPRAV NFGAAGSIMA HELLHIFYQL LLPGGCLACD NHALQEAHLC LKRHYAAFPL PSRTSFNDL TFLENAADV G LAIALQAYS KRLLRHHGET VLPSLDLSPQ QIFFRSYAQV MCRKPSPQDS HDTHSPPHLR VHGPLSSTPA FARYFRCARG ALLNPSSRCQ

image



Product Name: Recombinant Human KEL Protein
Catalog #: PEH6677



Background

Alternative Names

Kell blood group glycoprotein (EC 3.4.24.-) (CD antigen CD238)

Note

For Research Use Only , Not for Diagnostic Use