


Product Name: Recombinant Human BAK1 Protein
Catalog #: PEH6682



Summary

Name	BAK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q16611
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQECGEPALP SASEEQVAQD TEEVFRSYVF YRHQQEQEAE GVAAPADPEM VTLP LQPSST MGQVGRQLAI IGDDINRRYD SEFQTMLQHL QPTAENAYEY FTKIATSLFE SGINWGRVVA LLGFGYRLAL HVYQHGLTGF LGQVTRFVVD FMLHHCIARW IAQRGGWVAA

image



Background

Alternative Names Bcl-2 homologous antagonist/killer (Apoptosis regulator BAK) (Bcl-2-like protein 7) (Bcl2-L-7)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BAK1 Protein
Catalog #: PEH6682

