


Product Name: Recombinant Human POLM Protein
Catalog #: PEH6684



Summary

Name	DPOLM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NP87
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EILAEAAGFE GSEGRLLTFC RAASVLKALP SPVTTLSQLQ GLPHFGEHSS RVVQELLEHG VCEEVERVRR SEREQTMKLF TQIFGVGVKT ADRWYREGLR TLDDLREQPQ KLTQQQKAGL QHHQDLSTPV LRSDVDALQQ VVEEAVGQAL PGATVTLTGG FRRGKLQGHG VDFLITHPKE

image



Background

Alternative Names	DNA-directed DNA/RNA polymerase mu (Pol Mu) (EC 2.7.7.7) (Terminal transferase)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human POLM Protein
Catalog #: PEH6684

