


Product Name: Recombinant Human DUSP11 Protein
Catalog #: PEH6691



Summary

Name	DUS11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75319
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSSDFEDSAH LMQPVHNKPV KQGPRYNLHQ IQGHSAPRHF HTQTQSLQQS VRKFSENPBV YQRHHLPPPG PPGEDYSHRR

image



Background

Alternative Names	RNA/RNP complex-1-interacting phosphatase (EC 3.1.3.-) (Dual specificity protein phosphatase 11) (Phosphatase that interacts with RNA/RNP complex 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use