


Product Name: Recombinant Human ABCB5 Protein
Catalog #: PEH6698



Summary

Name	ABCB5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q2M3G0
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVIDKKPSID NFSTAGYKPE SIEGTVEFKN VSFNYPSPRS IKILKGLNLR IKSGETVALV GLNGSGKSTV VQLLQRLYDP DDGFIMVDEN DIRALNVRHY RDHIGVVSQE PVLFGTTISN NIKYGRDDVT DEEMERAARE ANAYDFIMEF PNKFNTLVGE KGAQMSGGQK QRIAIARALV RNPKILILDE ATSALDSESK SAVQAALEKA SKGRTTIVVA HRLSTIRSAD LIVTLKDGML AEKGHAELM AKRGLYSLV MSQDIKKADE QMESMTYSTE RKTNSLPLHS VKSIKSDFID KAEESTQSKE ISLPEVSLLK ILKLNKPEWP

image



Background

Alternative Names	ATP-binding cassette sub-family B member 5 (ABCB5 P-gp) (P-glycoprotein ABCB5)
--------------------------	--

Product Name: Recombinant Human ABCB5 Protein
Catalog #: PEH6698



Note

For Research Use Only , Not for Diagnostic Use