


Product Name: Recombinant Human ADCY2 Protein
Catalog #: PEH6705



Summary

Name	ADCY2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08462
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EREEIETMEN LNRVLLENVL PAHVAEHFLA RSLKNEELYH QSYDCVCVMF ASIPDFKEFY TESDVNKEGL ECLRLLNEII ADFDILLSKP KFSGVVEKIKT IGSTYMAATG LSAVPSQEHS QEPERQYMHI GTMVEFAFAL VGKLDAINKH SFNDFKLRVG INHGPIAGV IGAQKPQYDI WGNVTNVASR MDSTGVLDKI QVTEETSLVL QTLGYTCTCR GIINVKGKGD

image



Background

Alternative Names Adenylate cyclase type 2 (EC 4.6.1.1) (ATP pyrophosphate-lyase 2) (Adenylate cyclase type II) (Adenylyl cyclase 2)

Note

Product Name: Recombinant Human ADCY2 Protein
Catalog #: PEH6705



For Research Use Only , Not for Diagnostic Use