


Product Name: Recombinant Human COL1A2 Protein
Catalog #: PEH6714



Summary

Name	CO1A2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08123
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGSKGESGNK GEPGSAGPQG PPGPSGEEGK RGPNGEAGSA GPPGPPGLRG SPGSRGLPGA DGRAGVMGPP GSRGASGPAG VRGPNGDAGR PGEPGLMGPR GLPGSPGNIG PAGKEGPVGL PGIDGRPGPI GPAGARGEPP NIGFPGPKGP TGDPGKNGDK

image



Background

Alternative Names Collagen alpha-2(I) chain (Alpha-2 type I collagen)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human COL1A2 Protein
Catalog #: PEH6714

