

Product Name: Recombinant Human CHST1 Protein
Catalog #: PEH6717



Summary

Name	CHST1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O43916
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ADRRVMLGAS RDLLRSLYDC DLYFLENYIK PPPVNHTTDR IFRRGASRVL CSRPVCDPPG PADLVLEEGD CVRKCGLLNL TVAAEACRER SHVAIKTVRV PEVNDLRALV EDPRLNLKVI QLVRDPRGIL ASRSETFRDT YRLWRLWYGT GRKPYNLDVT QLTTVCEDFS NSVSTGLMRP PWLKGKGYMLV RYEDLARNPM

image



Background

Alternative Names	Carbohydrate sulfotransferase 1 (EC 2.8.2.21) (Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 1) (GST-1) (Keratan sulfate Gal-6 sulfotransferase) (KS6ST) (KSGal6ST) (KSST)
--------------------------	---

Note

Product Name: Recombinant Human CHST1 Protein
Catalog #: PEH6717



For Research Use Only , Not for Diagnostic Use