


Summary

Name	CHST4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8NCG5
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSLFQWENSR ALCSAPACDI IPQDEIIPRA HCRLLCSQQP FEVVEKACRS YSHVVLKEVR FFNLQSLYPL LKDPSLNLHI VHLVRDPRAV FRSRERTKGD LMIDSRIVMG QHEQKLLKED

image



Background

Alternative Names	Carbohydrate sulfotransferase 4 (EC 2.8.2.-) (Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 3) (GST-3) (High endothelial cells N-acetylglucosamine 6-O-sulfotransferase) (HEC-GlcNAc6ST) (L-selectin ligand sulfotransferase) (LSST) (N-acetylglucosamine 6-O-sulfotransferase 2) (GlcNAc6ST-2) (Gn6st-2)
--------------------------	--

Product Name: Recombinant Human CHST4 Protein
Catalog #: PEH6719



Note

For Research Use Only , Not for Diagnostic Use