


Product Name: Recombinant Human COL6A5 Protein
Catalog #: PEH6724



Summary

Name	CO6A5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A8TX70
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDMKATCVNM TSPNPENGGT ENTVLLPGL YEIKTENGDL FDEFDSQAQH LLVLGNHSS GSETATDLMQ KLYLLFSTEK LAMKDKEKAH LEEISALVVD

image



Background

Alternative Names	Collagen alpha-5(VI) chain (Collagen alpha-1(XXIX) chain) (von Willebrand factor A domain-containing protein 4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use