



Product Name: Recombinant Human RPP38-DT Protein
Catalog #: PEH6728



Summary

Name	CJ111_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N326
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ENQDKREKEY GVKHMPMGNN AGNLEPEKRR AVRVALSSAT AAQNIPSSVH CGCSKQWRLR LPSESLQSRG

image



Background

Alternative Names Putative uncharacterized protein RPP38-DT (RPP38 divergent transcript protein)

Note

For Research Use Only , Not for Diagnostic Use