


Product Name: Recombinant Human AZU1 Protein
Catalog #: PEH6735



Summary

Name	CAP7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P20160
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IVGGRKARPR QFPFLASIQN QGRHFCCGAL IHARFVMTAA SCFQSQNPGV STVVLGAYDL RRRERQSRQT FSISSMSENG YDPQQNLNDL MLLQLDREAN LTSSVTILPL PLQNATVEAG TRCQVAGWGS QRSGGRLSRF PRFVNVTVTP EDQCRPNVC TGVLTRRGGI CNGDGGTPLV CEGLAHGVAS FSLGPCGRGP DFFTRVALFR DWIDGVLNNP GPGP

image



Background

Alternative Names Azurocidin (Cationic antimicrobial protein CAP37) (Heparin-binding protein) (HBP) (hHBP)

Note

Product Name: Recombinant Human AZU1 Protein
Catalog #: PEH6735



For Research Use Only , Not for Diagnostic Use