


Product Name: Recombinant Human EBF2 Protein
Catalog #: PEH6737



Summary

Name	COE2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HAK2
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKTNNGTHYK LQLLYSNGVR TEQDLYVRLI DSVTKQPIAY EGQKNPEMC RVLLTHEVMC SRCCEKKSCG

image



Background

Alternative Names Transcription factor COE2 (Early B-cell factor 2) (EBF-2)

Note

For Research Use Only , Not for Diagnostic Use