


Product Name: Recombinant Human GLCE Protein
Catalog #: PEH6741



Summary

Name	GLCE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94923
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CSDDKAIQFP RRSSSGFRVD GF EKRAAASE SNNYMNHVAK QQSEEAFPQE QQKAPPVVGG FNSNVGSKVL GLKYEEIDCL INDEHTIKGR

image



Background

Alternative Names	D-glucuronyl C5-epimerase (EC 5.1.3.17) (Heparan sulfate C5-epimerase) (Hsepi) (Heparin/heparan sulfate:glucuronic acid C5-epimerase) (Heparosan-N-sulfate-glucuronate 5-epimerase)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use