


Product Name: Recombinant Human BMP10 Protein
Catalog #: PEH6747



Summary

Name	BMP10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	O95393
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NAKGN YCKRT PLYIDFKEIG WDSWIIAPPG YEAYECRGVC NYPLAEHLTP TKHAI IQALV H LKNSQKASK ACCVPTKLEP ISILYLDKGV VTYKFKYEGM AVSECGCR

image



Background

Alternative Names Bone morphogenetic protein 10 (BMP-10)

Note

For Research Use Only , Not for Diagnostic Use